




MARGARET KISS, PhD
Director of R&D
mkiss@ancera.com
Ancera



New data on *Salmonella* growth rates question current testing methods

“Biases such as bacterial-growth kinetics, variations in selective enrichment media and response to antimicrobial interventions widely differ across *Salmonella* serotypes and strains and can impact our awareness of food-safety threats.”

Conventional diagnostic methods for detecting *Salmonella* in live poultry production typically rely on analyzing isolated bacterial colonies from enriched samples.

Recent evidence, however, suggests these methods can overlook the presence of certain strains and serotypes, underestimating the true prevalence and risk of *Salmonella* in a sample, says Margaret Kiss, PhD, director of R&D at Ancera.

To effectively reduce *Salmonella* risk in a poultry-production system, she says, companies need accurate intelligence on:

1. **Detection:** Is *Salmonella* present?
2. **Quantity:** How much *Salmonella* is present?
3. **Strains and serovars:** What *Salmonella* strains and serovars are present?
4. **Susceptibility:** What kills or inhibits it?

“Biases such as bacterial-growth kinetics, variations in selective enrichment media and response to antimicrobial interventions widely differ across *Salmonella* serotypes and strains and can impact our awareness of food-safety threats,” Kiss explains.

DIFFERENCES IN GROWTH KINETICS

Bacterial-growth kinetics are often described as having four phases: the lag phase, exponential phase, stationary phase and death phase.


Each is characterized by distinct growth rates and physiological states. Different *Salmonella* strains (genetic variants) and serovars (distinct variations within a species) exhibit different growth kinetics, which can impact the detectability of *Salmonella* in enriched samples.

“Unique characteristics in the length of the phases will impact the number of cells present after a given enrichment time,” Kiss adds.

continued

“ Unique characteristics in the length of the phases will impact the number of cells present after a given enrichment time. ”

Figure 1. Comparison of unstructured kinetic bacterial-growth models



Source: Mulaiwa M, Nyende-Byakika S, Dinka M. 2020

Figure 2. Lag time in strains of different serovars





Figure 3. Lag time in different strains of the same serovar



* Enrichment media = buffered peptone water (BPW) • + Starting inoculum between 10-100 CFU/ml



SEROTYPE AND STRAIN-SPECIFIC DIFFERENCES IN LAG TIME

Different *Salmonella* serovars and strains exhibit considerable variability in duration of the lag phase.

Ancera’s internal research shows that lag time varied from 6 to 9 hours at low starting inoculum (10-100 cfu/ml) in buffered peptone water (BPW). Even strains of the same serotype exhibited variability in lag times.

“Food-safety experts must incorporate differences in strain-specific lag time to properly interpret the results of *Salmonella* isolation and quantification,” Kiss says.

DOUBLING-TIME VARIABILITY

Doubling time — the amount of time for a bacterial population to double in size — also has significant variability among different *Salmonella* strains and serovars.

In Ancera’s research, doubling time varied between 19 and 30 minutes for the different strains and serovars examined in BPW.

“And in addition to differences in lag times, different strains of the same serotype had different doubling times,” Kiss reports. “For this set of strains and growth conditions, a quantitative *Salmonella* measurement could be impacted in a mixed serovar population.”

“Food-safety experts must incorporate differences in strain-specific lag time to properly interpret the results of *Salmonella* isolation and quantification.”

continued

Figure 4. Doubling time in strains of different serovars





Figure 5. Doubling time in different strains of the same serovar



* Enrichment media = buffered peptone water (BPW) • + Starting inoculum between 10-100 CFU/ml

“...the culture media can influence the *Salmonella*-growth kinetics, depending on the composition of the media and the characteristics of the *Salmonella* strains present in the sample.”

For example, after 9 hours the *Salmonella* Kentucky population would have gone through about six doublings while the *Salmonella* Infantis population would, on average, still be in lag phase. Accordingly, *S. Kentucky* would be over-represented in the sample relative to *S. Infantis* after enrichment.


SELECTIVE MEDIA IN SALMONELLA SURVEILLANCE

Selective media facilitates the isolation and identification of target organisms from samples. To help with detection, these media typically contain ingredients that inhibit the growth of non-target bacteria while allowing the growth of *Salmonella*.

“However, the culture media can influence the *Salmonella*-growth kinetics, depending on the composition of the media and the characteristics of the *Salmonella* strains present in the sample,” Kiss cautions.

Larsen et al. (2021) used CRISPR SeroSeq to demonstrate that the choice of culture media can influence the competitive dynamics between different strains. Figure 6 shows specific *Salmonella* Montevideo (white) and *Salmonella* Typhimurium (black) strains that were mixed at equal ratios and cultured in two different broths. In Rappaport-Vassiliadis soya peptone broth, the *S. Montevideo* strain outgrew the *S. Typhimurium* strain, whereas in tetrathionate broth, the *S. Typhimurium* strain outgrew the *S. Montevideo* strain.

Figure 6. Mixed *Salmonella* cultures reveal competitive advantages between strains during pre-enrichment and selective enrichment.



Larsen BR, Richardson KE, Obe T, Schaeffer C, Shariat NW. J Food Safety. 2021;41(6). e12934.

ANTIMICROBIAL SUSCEPTIBILITY AMONG SALMONELLA STRAINS

Antimicrobial susceptibility is another important consideration in implementing a *Salmonella*-intervention process in live production.


Ancera’s research showed similar growth kinetics among three different *Salmonella* strains in rich media in the absence of any antimicrobial.

However, in the presence of 312 ppm peracetic acid (PAA) — an organic chemical used to control *Salmonella* in processing — the results revealed different responses among the strains, with some strains resuming growth shortly after exposure to PAA, while others remained unaffected even after extended incubation periods, Kiss says.

A strain of *S. Infantis* resumed growing within 4 hours of exposure, while a strain



Figure 7. Antimicrobial susceptibility of *S. Enteritidis*, *S. Infantis* and *S. Kentucky*



of *S. Enteritidis* resumed growing within 8 hours after exposure and a strain of *S. Kentucky* did not resume growth at all within 12 hours.

SALMONELLA MONITORING AT SCALE

Proper interpretation of *Salmonella*-surveillance results requires a deep understanding of growth kinetics, selective media and antimicrobial susceptibility. According to Kiss, there is inherent complexity within the variability in lag time, doubling time and response to selective media among different strains and serovars.

As food-safety professionals develop multi-layered food-defense systems that effectively detect, quantify, categorize and intervene against *Salmonella*, it is

critical that they have granular serotype information and robust data systems to understand and explain the interconnected relationship between varying strains.

“Food producers looking for a more detailed understanding of their *Salmonella*-control operations will require interdisciplinary collaborations between microbiologists, epidemiologists and food-safety experts,” Kiss says.

“They will need to leverage advanced molecular biology, genomics and data analytics. They will need researchers who can further improve the sensitivity, specificity and efficiency of *Salmonella*-surveillance methods. Or they can work with Ancera, which does all of this and more.”

continued

New data on *Salmonella* growth rates question current testing methods

“ Food producers looking for a more detailed understanding of their *Salmonella*-control operations will require interdisciplinary collaborations between microbiologists, epidemiologists and food-safety experts. ”

For more information, contact Margaret Kiss, PhD, at mkiss@ancera.com or visit ancera.com.

